

## Instructions for use

### **Attention!!**

**Arteriosclerosis, or who are suspected, please do not use this product.**  
**It may take an extra burden on the arterial vessel.**

#### 1. Adjusting length of the belt

**Note:** Please do not cut the belt if you are planning to return the belt within the trial period.

Please cut off the belt with scissors so that the end of the belt protrudes at least 8 cm from the buckle. ( $A \geq 8$  cm)

When wearing HAIR REGROWTH BELT to your head as shown in the lower right figure 1. It is difficult to tighten the belt strongly if it is too short.

#### 2. Wearing process


Put it around your head as shown in Fig. 1 finally by the following wearing process.

(You can see a [video](#) on the web site of HAIR REGROWTH BELT.


URL [http://www.ip-enterprisez.com/way\\_of\\_wearing\\_E](http://www.ip-enterprisez.com/way_of_wearing_E)

- ① First, put the belt around your head in a loose state. If your forelock is long, fasten them with a hair clip etc. so that the forelock does not hang over the forehead, or your forelock will be caught in the buckle.  
 ※ 《Important》 Make sure that the smaller arc of the belt is the top.
- ② Pull out the hairs whose roots are located on the upper side (the vertex side) than the wearing belt from the back of the belt using a rattail comb. As a result of that, the upper inside of the belt contacts with your scalp and the belt is prevented to slip and slide up toward your vertex.
- ③ As shown in Fig. 2, hold the end of the belt with the right hand, push the tip of the right thumb against the buckle and pull the right hand to the right rearward (arrow direction) using it as a fulcrum, you can tighten it strongly. Please close the lid of the buckle with your left hand in that state. When shutting the lid of buckle, be careful so that your finger wouldn't be caught up. Please close the lid of the buckle surely.
- ④ The wearing position of the belt shall be such that the belt naturally does not rise after wearing and the tightening force does not become loose.

**Fig. 1**



**Fig. 2**



### Criterion of tightening force

**When going to bed after wearing :**  
**strongly (enough to be able to wrinkle up the scalp of the crown )**

**In other cases :**  
**weakly (The extent to which the stretched crown scalp loosens )**

**It is basic to make as high a tightening force as possible within the range not feeling the pain due to the tightening of the belt. According to this basic, if you feel pain, you can adjust the tightening force by loosening the tightening force. To slightly loosen the tightening force, please shift the belt on the forehead upward.**

### 3. About the frequency of wearing and the time of wearing

Wearing the belt every day is recommended. At least more than 3 hours are needed at once of wearing time. Therefore, going to bed (lying down) while wearing the belt is recommended.

When waking up in the next morning, please take off the belt. If the belt slips and loosens when waking up, next time please adjust the wearing position of the belt lower slightly than the previous time.

When sleeping in bed after wearing the belt, lie down as soon as possible (within 15 minutes) after wearing. When more than 30 minutes are passed without lying down, vertex scalp falls in the congestion, and you may feel the pain by congestion, so please be careful.

### 4. To maximize the effect of HAIR REGROWTH BELT

Wearing the belt "long time as much as possible with tightening force as high as possible" is recommended. Therefore, going to bed (lying down) while wearing the belt is recommended. However, please control the tightening force to not feel pain due to tightening of the belt.

### 5. Notes on first use

Wearing the belt leaves a mark on your forehead. It takes about one hour for the mark not to stand out. Wearing at the night before a holiday when you have no appointment to see others and go out, make sure how long it takes for the mark not to stand out. And then, control the time to see others and go out. Further, you may have bed hair by wearing the belt, so time to fix the bed hair is also needed.

### 6. Pulsating of scalp

Pulsating around the vertex may be felt while wearing the belt and also without wearing. It is normal and a proof that the blood circulation of the scalp is improved by the effect of the belt. There is no need to worry.

### 7. In case of feeling a pain

If you feel pain due to tightening of the belt during use, loosen the tightening force of the belt or remove the belt from the head. When using it again, please try to weaken the tightening force a little so that you can use it without feeling pain.

n.

### 8. About allergy.

Occasionally, some people have an allergy to metal which is made of the buckle. Please stop using the belt in this case and consult a doctor. Hypoallergenic metal buckle is available for metal allergy. Parts of the buckle in contact with a forehead are coated with thin layer (pencil hardness 5H or more) of silica and titanium oxide harmless to human body. If you worry about metal allergy, please buy this.

### 9. About the durability limit

Because the belt is made of genuine leather, it gradually extends when you use it. The effect will diminish as the belt stretches and the curvature of the belt decreases. If you feel that the effect has diminished and hair loss has increased, please purchase a new belt. It is also possible to purchase only replacement belts.

### 10. Maintenance

The belt is made of genuine leather. Please do not wash it or wet it with water. There is no problem for the belt to sweat. Please dry away from direct sunlight everyday. When washing the buckle, please brush with a flexible brush. Please do not brush strongly a buckle for allergy as the coated layer may be peeled off.

### 11. Others

Please do not use it when you are in bad physical condition. Even if the tightening force is weak, it becomes easy to feel pain when you are in poor physical condition.

Since the belt is made of genuine leather, there are wrinkles, spots, small scratches, etc., due to the difference in skin part. Since there is no influence on the function / effect of the belt, please use it as it is.

When disposing, please follow the rules of the local dump.

### 12. Specification of HAIR REGROWTH BELT

Belt (made in Japan): Made of cow leather (Raw material may be foreign-made, even in that case the tanning is in Japan.)

Buckle (Made in Japan): Stainless steel (SUS 304)

(Anti-allergy layer: thin layer of titanium oxide and silica)

Patented in Japan (No. 5255440).

#### 《Separate sale items》

Replacement belt: Please use for replacing the stretched belt.

Buckle for metal allergy : Parts of the buckle in contact with a forehead are coated with thin film (pencil hardness 5H or more) of silica and titanium oxide harmless to human body.